

# KeyTec® TR-FRET

## mAb anti-HA-Solar Eu



CAT. & Size    A1020017S (1,000 tests)  
                    A1020017L (10,000 tests)  
Storage at     -60°C or below

VKEYBIO-01-2024  
For Research Use Only  
Not For Diagnostic Or Therapeutic Use

## KeyTec® TR-FRET

### mAb anti-HA-Solar Eu

#### Instruction Manual

#### 1. Introduction

KeyTec® TR-FRET mAb anti-HA-Solar Eu is designed for developing the TR-FRET Assay. The anti-HA antibody is a mouse monoclonal antibody. In the Protein-Protein Interaction assay, one HA-tagged protein binds to the donor (KeyTec® TR-FRET mAb anti-HA-Solar Eu<sup>\*1</sup>), and the other protein is labeled (directly or indirectly) with the acceptor (KeyTec® TR-FRET LA/HX<sup>\*2</sup>). When the two proteins interact, the donor molecule is brought into proximity with the acceptor molecule. Excitation of the donor will result in the generation of the TR-FRET signal at 665 nm, proportional to the extent of protein interaction.

\*1 KeyTec® TR-FRET Solar Eu: TR-FRET Donor Molecule

\*2 KeyTec® TR-FRET LA/HX: TR-FRET Acceptor Molecule




Figure 1. KeyTec® TR-FRET Protein-Protein Interaction assay mode

## 2. Components

| Components                  | A1020017S<br>(1,000 tests) | A1020017L<br>(10,000 tests) |
|-----------------------------|----------------------------|-----------------------------|
| KeyTec® TR-FRET             | 1 vial                     | 1 vial                      |
| mAb anti-HA-Solar Eu (100X) | 50 µL/vial                 | 500 µL/vial                 |

| KeyTec® Materials Required But Not Supplied   | CAT. & Size                |
|---|----------------------------|
| KeyTec® TR-FRET Binding Assay Diluent Buffer  | A1010001L<br>(200 mL)      |
| KeyTec® TR-FRET Solar Eu Detection Buffer   | A1010002L<br>(120 mL)      |
| KeyTec® 384-Well White Flat Low-Volume Microplates,<br>PS, Solid, Non-treated, No lid | M2000102N<br>(40 Pcs/Box)  |
| KeyTec® Fluorescent High-Transparency Microplate Top Seals                            | M1000102N<br>(100 Pcs/Box) |

## 3. Storage Conditions

- Upon receipt, store the reagent below -60 °C
- Up to 1 years from date of receipt, when stored and handled as recommended.
- When first thaw, aliquot the reagents as needed to avoid multiple freeze-thaw cycles.

## 4. Assay Procedure

### 1.1 Assay Format

| Assay Format  | Total Volume (20 µL <sup>3</sup> ) |
|---|------------------------------------|
| Other assay components                                    | 10 µL                              |
| KeyTec® TR-FRET Donor (Solar Eu/Tb) working solution (1X) | 5 µL                               |
| KeyTec® TR-FRET Acceptor (LA/HX) working solution (1X)    | 5 µL                               |

\*<sup>3</sup> The assay volume is optimized for 384-well microplates, and can be adjusted proportionally to perform in 96- or 1536-well microplates.

## 1.2 Reagents Handling

### 1) Buffers

- ◆ KeyTec® TR-FRET Solar Eu Detection Buffer (A1010002L) has been optimized for maximum performance.
- ◆ Use the same buffer to prepare both the donor and the acceptor (LA/HX) conjugates.
- ◆ KeyTec® TR-FRET Binding Assay Diluent Buffer (A1010001L) is recommended for dilution and preparation of other components or samples.
- ◆ If using a homemade buffer solution, avoid SDS and ensure KF addition.

### 2) Conjugates

- ◆ Thaw reagents on ice and equilibrate to room temperature before use.
- ◆ Prepare working solutions as per the purchased product instructions. The storage solution for KeyTec® TR-FRET mAb anti-HA-Solar Eu is 100X; dilute 100 times for a 1X working solution. For example, mix 50 µL of the storage solution with 4950 µL of KeyTec® TR-FRET Solar Eu Detection Buffer for a 1X working solution.
- ◆ Optimal amounts per well can be further optimized based on different assay format and conditions.

## 1.3 Data Calculating

- ◆ Calculate the ratio of 665 nm/615 nm (TR-FRET Ratio) and the CV for each individual well.

$$\text{TR-FRET Ratio} = \frac{\text{Signal 665 nm}}{\text{Signal 615 nm}} \times 10,000$$