

KeyTec® TR-FRET

pAb anti-Mouse IgG-Solar Eu



CAT. & Size A1020029S (1,000 tests)
 A1020029L (10,000 tests)
Storage at -60°C or below

VKEYBIO-01-2024
For Research Use Only
Not For Diagnostic Or Therapeutic Use

KeyTec® TR-FRET

pAb anti-Mouse IgG-Solar Eu

Instruction Manual

1. Introduction

KeyTec® TR-FRET pAb anti-Mouse IgG-Solar Eu is designed for developing the TR-FRET Assay. The anti-Mouse IgG antibody is a goat polyclonal antibody. In the Antigen-Antibody Interaction assay, the Mouse IgG antibody binds to the donor (KeyTec® TR-FRET pAb anti-Mouse IgG-Solar Eu^{*1}), and the antigen is labeled (directly or indirectly) with the acceptor (KeyTec® TR-FRET LA/HX^{*2}). When the Antigen-Antibody interact, the donor molecule is brought into proximity with the acceptor molecule. Excitation of the donor will result in the generation of the TR-FRET signal at 665 nm, proportional to the extent of antigen-antibody interaction.

*¹ KeyTec® TR-FRET Solar Eu: TR-FRET Donor Molecule

*² KeyTec® TR-FRET LA/HX: TR-FRET Acceptor Molecule




Figure 1. KeyTec® TR-FRET Antigen-Antibody Interaction assay mode

2. Components

| Components | A1020029S (1,000 tests) | A1020029L (10,000 tests) |
|------------------------------------|----------------------------|-----------------------------|
| KeyTec® TR-FRET | 1 vial | 1 vial |
| pAb anti-Mouse IgG-Solar Eu (100X) | 50 µL/vial | 500 µL/vial |

| KeyTec® Materials Required But Not Supplied | CAT. & Size |
|---|----------------------------|
| KeyTec® TR-FRET Binding Assay Diluent Buffer | A1010001L (200 mL) |
| KeyTec® TR-FRET Solar Eu Detection Buffer | A1010002L (120 mL) |
| KeyTec® 384-Well White Flat Low-Volume Microplates, PS, Solid, Non-treated, No lid | M2000102N (40 Pcs/Box) |
| KeyTec® Fluorescent High-Transparency Microplate Top Seals | M1000102N (100 Pcs/Box) |

3. Storage Conditions

- Upon receipt, store the reagent below -60 °C
- Up to 1 years from date of receipt, when stored and handled as recommended.
- When first thaw, aliquot the reagents as needed to avoid multiple freeze-thaw cycles.

4. Assay Procedure

1.1 Assay Format

| Assay Format | Total Volume (20 µL ³) |
|---|------------------------------------|
| Other assay components | 10 µL |
| KeyTec® TR-FRET Donor (Solar Eu/Tb) working solution (1X) | 5 µL |
| KeyTec® TR-FRET Acceptor (LA/HX) working solution (1X) | 5 µL |

*³ The assay volume is optimized for 384-well microplates, and can be adjusted proportionally to perform in 96- or 1536-well microplates.

1.2 Reagents Handling

1) Buffers

- ◆ KeyTec® TR-FRET Solar Eu Detection Buffer (A1010002L) has been optimized for maximum performance.
- ◆ Use the same buffer to prepare both the donor and the acceptor (LA/HX) conjugates.
- ◆ KeyTec® TR-FRET Binding Assay Diluent Buffer (A1010001L) is recommended for dilution and preparation of other components or samples.
- ◆ If using a homemade buffer solution, avoid SDS and ensure KF addition.

2) Conjugates

- ◆ Thaw reagents on ice and equilibrate to room temperature before use.
- ◆ Prepare working solutions as per the purchased product instructions. The storage solution for KeyTec® TR-FRET pAb anti-Mouse IgG-Solar Eu is 100X; dilute 100 times for a 1X working solution. For example, mix 50 µL of the storage solution with 4950 µL of KeyTec® TR-FRET Solar Eu Detection Buffer for a 1X working solution.
- ◆ Optimal amounts per well can be further optimized based on different assay format and conditions.

1.3 Data Calculating

- ◆ Calculate the ratio of 665 nm/615 nm (TR-FRET Ratio) and the CV for each individual well.

$$\text{TR-FRET Ratio} = \frac{\text{Signal 665 nm}}{\text{Signal 615 nm}} \times 10,000$$