

CAT. & Size M2000502N (40 Pcs/Box)

VKEYBIO-01-2024

Storage at RT

For Research Use Only

Not For Diagnostic Or Therapeutic Use

KeyTec® 1536-Well White Flat Microplates


PS, Solid, Non-treated, No lid

Instruction Manual

1. Introduction

KeyTec® 1536-Well White Flat Microplates is designed for TR-FRET, Luminescent, Alpha and radiometric assays. With shallow wells designed for lower assay volumes (suitable for 4 μ L–12 μ L), it brings samples in closer proximity to the reader's detectors. Designed to miniaturize your assay and utilize cost reduction advantages.

2. Dimension



3. Specifications

CAT.	Description
Automation Compatible	Yes
Coating Treatment	Untreated
Color	White
Detection Method	Luminescence, Alpha, Time-Resolved Fluorescence (TRF & TR-FRET), Radiometric
Material	Polystyrene
Size	40 Pcs/Box
Wells Number	1536 well plate
Well Shape	Shallow-well
Working Volume	4 μ L – 12 μ L
Total Volume	12 μ L