

CAT. & Size M2000602N (40 Pcs/Box)

VKEYBIO-01-2024

Storage at RT

For Research Use Only

Not For Diagnostic Or Therapeutic Use

KeyTec® 1536-Well Black Flat Microplates


PS, Solid, Non-treated, No lid

Instruction Manual

1. Introduction

KeyTec® 1536-Well Black Flat Microplates is designed for Fluorescence, Fluorescence Polarization, Luminescent assays. With shallow wells designed for lower assay volumes (suitable for 4 μ L–12 μ L), it brings samples in closer proximity to the reader's detectors. Designed to miniaturize your assay and utilize cost reduction advantages.

2. Dimension



3. Specifications

CAT.	Description
Automation Compatible	Yes
Coating Treatment	Untreated
Color	Black
Detection Method	Fluorescence, Fluorescence Polarization, Luminescent
Material	Polystyrene
Size	40 Pcs/Box
Wells Number	1536 well plate
Well Shape	Shallow-well
Working Volume	4 μ L – 12 μ L
Total Volume	12 μ L