

CAT. & Size M2000302N (40 Pcs/Box)

VKEYBIO-01-2024

Storage at RT

For Research Use Only

Not For Diagnostic Or Therapeutic Use

## KeyTec® 384-Well White Flat Microplates


### PS, Solid, Non-treated, No lid

#### Instruction Manual

#### 1. Introduction

KeyTec® 384-Well White Flat Microplate is designed for TR-FRET, Luminescent, Alpha and radiometric assays. With square wells designed for standard assay volumes (suitable for 20 - 100  $\mu$ L), it ensures easier operation and reduces errors.

#### 2. Dimension



### 3. Specifications

CAT.	Description
Automation Compatible	Yes
Coating Treatment	Untreated
Color	White
Detection Method	Luminescence, Alpha, Time-Resolved Fluorescence (TRF & TR-FRET), Radiometric
Material	Polystyrene
Size	40 Pcs/Box
Wells Number	384 well plate
Well Shape	Square well
Working Volume	20 $\mu$ L - 100 $\mu$ L
Total Volume	120 $\mu$ L